


Complete control with **SIMPLANT MASTER**

SimPlant Master allows for complete control in all aspects of SimPlant. Read CT and conebeam scanner data directly into SimPlant Master, create your own 3D images and convert files into all other SimPlant software types for your colleagues.



FEATURES INCLUDE:

- *Convert Master files into Pro, Planner and View*
- *Read data from any type of Conebeam or CT scanner*
- *Highlight the ID canal*
- *Virtually place real implants*

- *Access a comprehensive implant library*
- *Produce fully interactive 3D images in-house*
- *Link to computer guided surgery [SurgiGuides from Materialise]*



COMPUTER GUIDED
IMPLANTS LTD
3B Kitsbury Rd
Berkhamsted
Herts HP4 3EG

T 01442 872787
F 01442 877630

sales@cgimplants.co.uk
www.cgimplants.co.uk