


The Classic: SIMPLANT PLANNER

SimPlant Planner is a precise, pre-operative software, which allows for real 3D planning – assisting you in finding the perfect position for each implant – so you can achieve an accurate and aesthetically pleasing result.

FEATURES INCLUDE:

- Fully interactive 3D image
- Highlight the ID canal
- Virtually place real implants
- Link to computer guided surgery [SurgiGuides from Materialise]



- Access a comprehensive implant library
- Panoramic, Axial and Cross-sectional views
- Make 2D and 3D bone measurements
- And more!!



COMPUTER GUIDED
IMPLANTS LTD
3B Kitsbury Rd
Berkhamsted
Herts HP4 3EG
T 01442 872787
F 01442 877630
sales@cgimplants.co.uk
www.cgimplants.co.uk